
4/1990

**Gebrauchsanweisung
Instruction Sheet
Mode d'emploi**
348 18

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The large gyroscope is used to study all the fundamental laws of gyroscope action and to demonstrate the rotary motions by means of the gyroscope, e.g. the neutral equilibrium gyroscope.

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Literature:
New Physics Leaflets for Colleges and Universities, Vol. 1 (599 952)

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1. Description

This gyroscope is shaped like a bicycle wheel with spokes. It has a small weight, but a large moment of inertia. The hub has the shape of a bee-hive and is provided with a ball-bearing. This contraption leaves the centre of gravity of the gyroscope accessible, and the gyroscope can be supported above, at, and below its centre of gravity. Only a gyroscope in which the point of support can be varied in this manner will permit a demonstration of all the fundamental laws of gyroscopic action.

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A rod of 50 cm length is supplied with the gyroscope for performing experiments in which the gyroscope is held by hand. In place of this long rod, a short rod with pointed end can be introduced into the ball-bearing of the gyroscope and fastened by means of the lateral thumb screw. Another short rod with a cup on top can be clamped in the large tripod stand supplied with

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the apparatus and will then provide a support for the gyroscope.

The rod fastened in the ball-bearing may be displaced in the longitudinal direction by loosening the thumb screw. Thus experiments may be carried out at different positions of the point of support relative to the centre of gravity, and it is even possible to make such changes while the gyroscope is in rotation. When the gyroscope is supported at its centre of gravity, it will rotate about the pointed end of the short rod; in all other cases it will rotate about the ball-bearing, at almost equally small friction.

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2. Operation

The gyroscope is very easy to start: it should be firmly held in the left hand by a rod inserted into the ball-bearing, and given a brief, but large angular momentum by a forceful push against one of the spokes with the right. Suppose, for example, the experimenter wishes to arrange the spinning gyroscope on the rod with cup in the tripod stand. He should then insert the pointed rod, hold the gyroscope in the left by the part of the rod that protrudes above the ball-bearing, start the gyroscope with the right as described above, and set the point on the cup. This will not be found difficult after practising a few times.

Beyond this, the gyroscope can be operated by constant angular momentums at places of different diameters, while the axis is positioned horizontally and fastened by a wound-up cord (e.g. fishing line 309 48). This can be accomplished by suspending a few loads at the line.

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3. Experiments

Experiments in which the gyroscope is in neutral equilibrium

Fasten the pointed rod in the ball-bearing in a position where the annular mark on the rod is just above the upper edge of the ball-bearing. In that case, the gyroscope will be supported at its exact centre of gravity, and will therefore be in neutral equilibrium. This can be shown by placing the non-rotating gyroscope on the cup, making minor corrections if necessary. At this adjustment the axis of the spinning gyroscope will remain stationary in space if no additional angular momentum is acting on it, which might have been produced on adapting it on the cup, and the resultant torque on the gyroscope will be zero.

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If a gyroscope that rotates steadily is given an additional angular momentum by pushing against its "axis of symmetry" (that is the axis about which the moment of inertia of the gyroscope is greatest; it corresponds to the rod inserted into the ball-bearing in this example) it will perform an additional rotary motion called "nutation". In the case of this motion it is possible to show that the axis of symmetry and the "instantaneous axis of spin" (that is the axis about which the gyroscope happens to rotate at any given moment) both revolve about the "angular momentum axis" (that is the direction of the angular momentum vector). If printed newspaper is glued on the top side of the ball-bearing, the instantaneous axis of spin or the angular momentum vector can be made visible, respectively. At the point of intersection of the instantaneous axis of spin, individual letters will be seen, but everywhere else they will merge with the white of the paper to a uniform grey. The instantaneous axis of spin migrates in a circle around the angular momentum vector.

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Experiments in which an external force is acting on the gyroscope

Starting with the gyroscope in neutral equilibrium spinning without nutation it can be shown that any force acting at right angles to the axis of spin will cause this latter to move in the direction of the third coordinate. The easiest way of letting a permanent force of this kind act upon the gyroscope is to use the weight of the gyroscope as force and produce a moment arm by slightly displacing the point of support from the centre of gravity. To this end the upper thumb screw of the ball-bearing should be gripped while the gyroscope is spinning, and the gyroscope should be lifted or lowered somewhat to shift its point of support from the centre of gravity. A position is easy to find at which the gyroscope moves practically nutationless. It will then perform rotary motions on a cone-shaped shell, which are called precession.

The faster the gyroscope is spinning, the slower will be its precessional motion, the angular velocity of the precession being inversely proportional to the angular velocity of the rotary motion of the gyroscope. Moreover the angular velocity of the precession is proportional to the moment of the applied force, as can be shown by repeatedly varying the point of support of the gyroscope. The precessional direction will be

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inversed if the point of support of the gyroscope is moved from one side of the centre of gravity of the gyroscope to the other by shifting the pointed rod. The presumable direction of the precession can even be shown on a gyroscope supported at its centre of gravity if a rod is used to impart suitable angular momentums to it by hand.

A gyroscope not supported at its centre of gravity can perform nutational and precessional motions at the same time, so that very complex motions may result.

Other gyroscope experiments

If the 50 cm long rod is inserted into the gyroscope, other very important experiments can be conducted. Do not fail to hand the spinning gyroscope to students so that they may feel and experience for themselves the tendency of the gyroscope to "give way" in a direction normal to the applied force.

In conjunction with the rotatable disc (331 66) the theorem of the conservation of angular momentum can be verified in a very impressive manner. Standing or sitting (stool 331 69) on the rotating disc and holding the rapidly rotating gyroscope so that its axis of spin is vertical, then turning the gyroscope through 180°, an angular momentum will be transferred to the rotatable disc which will be twice as large as the one produced when stopping the gyroscope when its axis is vertical.

The large gyroscope can also be used to demonstrate the action of a gyrocompass. The gyroscope should be started in rapid rotation and then held by hand, as close to the centre of gravity as possible. The gyroscope will then tend to bring its axis of spin in a vertical position if oneself is turning as evenly as possible. This experiment is also preferably performed with the rotatable disc (331 66). In this model the rotating stool would correspond to the rotating earth.

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