
2/97-Di/Sf-

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The motion transducer (337 63), which consists of the motion sensing element (337 631) and the measurement transducer, converts the path  $s$ , the velocity  $v$  and the acceleration  $a$  of a body in motion into proportional voltages. This enables the instantaneous values of the characteristic quantities of a motion to be displayed or recorded with the correct sign using analog voltage meters (pointer instruments, recorders).

The quantities  $s$ ,  $v$  and  $a$  can be investigated both as a function of the time  $t$  and in terms of their relationship to one another. A built-in time base makes possible the use of XY recorders for registering time-dependent functions as well.

For computer-assisted experimenting, the digital information regarding the course of the motion supplied by the motion sensing element (337 631) can be input into the computer via a CASSY-interface (e.g. from 524 007) equipped with the BMW-box (524 032) – without conversion to analog voltages. Using specially optimized software, the information is processed so that the quantities path, time, velocity and acceleration are displayed correctly in terms of numerical value and dimension.

### Measuring principle

The motion of a body is transmitted to the motion sensing element (337 631) via an extremely low-friction spoked wheel, which is scanned by two miniature light barriers. The obscurations of the light barriers are counted by the measuring electronics of the measurement transducer or by the computer as the measure of the motion and converted into a path-proportional voltage or (after processing of the digital signals in the computer) displayed as the path. The velocity is determined from the number of obscurations per unit of time, and the acceleration from the change in this counting rate over time.

### Experiment examples:

- Linear motions, e.g. on the linear air track (337 50)
- Rotational motions, with the rotation model (347 23)
- Mechanical oscillations up to approx. 1 Hz, e.g. linear oscillations of spring and string pendulums, harmonic and chaotic rotational oscillations (e.g. with the torsion pendulum 346 00)

### Literature:

"Motion Transducer" (337 641, in German)  
"Measuring and Evaluating", physics part 2 (524 131, in German)  
Manual for CASSY software "Motion, BMW" (524 702)

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## 1 Safety notes

- Convert the measurement transducer as described in section 4.2 before using with a supply voltage of  $115 \text{ V AC}$ .
- Always unplug the device before opening the housing.
- Do not place loads greater than 50 g on the spoked wheel ① (greater loads cause excessive bearing friction and premature wear of the needle bearing).
- Do not apply any external voltage to outputs ⑧ - ⑪.
- Maximum permissible load of interrupter socket pair ⑬;  $24 \text{ V/1 A}$ .

## 2 Scope of supply, description, technical data

- ① - ④ Motion sensing element (337 631)
- ① Spoked wheel as motion sensor; mounted on low-friction needle bearings between two miniature light barriers; effective diameter (thread groove): 29.3 mm  
Equivalent mass of pulley: 1 g
- ② Socket for multicore cable ⑮ (501 16) for connecting transducer or (when processing the digital information supplied by the motion sensing element in a computer) to the BMW-box (524 032)
- ③ Threaded holes for screwing in stand rod ④
- ④ Stand rod (12 cm x 1 cm dia.), with thread for ③

Housing dimensions: 8 cm x 3 cm dia.

Weight: 0.15 kg

- ⑤ - ⑬ Measurement transducer with measuring electronics, voltage supply for the light barriers of the motion sensing element (337 631) and time base for controlling an XY-recorder
- ⑤ On-off switch with indicator lamp
- ⑥ Socket for multicore cable ⑮ (501 16) for connecting motion sensing element (337 631)
- ⑦ Reset key for s-output ⑧
- ⑧ - ⑪ Analog outputs (short-circuit proof) for connecting voltmeters, recorders or the analog inputs of an interface device (computer measuring station)  
Bottom sockets connected to each other and protective earth;  
operate delay < 0.2 s at ⑧ - ⑩; undelayed at ⑪  
Error:  
maximum:  $\pm 2\%$   
typical:  $\pm 1\%$
- ⑧ s-output with path-proportional voltage;  $1 \text{ V} \triangleq 1 \text{ m}$   
Measuring range from  $-1.5 \text{ m}$  to  $+2.5 \text{ m} \triangleq -1.5 \text{ V}$  to  $+2.5 \text{ V}$   
Internal resistance:  $100 \Omega$   
Resolution: better than 2 mm
- ⑨ v-output with velocity-proportional voltage;  $1 \text{ V} \triangleq 1 \text{ m s}^{-1}$   
Measuring range up to  $\pm 5 \text{ m s}^{-1} \triangleq \pm 5 \text{ V}$   
Internal resistance:  $100 \Omega$
- ⑩ a-output with acceleration-proportional voltage;  $1 \text{ V} \triangleq 1 \text{ ms}^{-2}$   
Measuring range up to  $\pm 10 \text{ m s}^{-2} \triangleq \pm 10 \text{ V}$   
Internal resistance:  $100 \Omega$
- ⑪ t-output with time-proportional voltage (time base)  $0.1 \text{ V} \triangleq 1 \text{ s}$   
Measuring range to  $120 \text{ s} \triangleq 12 \text{ V}$   
Internal resistance:  $1 \text{ k}\Omega \pm 5\%$
- ⑫ Control switch for time base ⑪ and interrupter sockets ⑬  
Start: voltage at socket pair ⑪ increases by  $0.1 \text{ V/s}$   
Stop: voltage constant  
Reset: voltage drops to  $0 \text{ V}$
- ⑬ Interrupter sockets, coupled with control switch ⑫; for interrupting circuits (e.g. with holding magnet) synchronous with start of time base.  
Contact between the sockets closed at reset, open at stop and start.  
Maximum permissible load:  $24 \text{ V/1 A}$ .
- ⑭ Perlon thread

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
- ⑮ Multicore cable, 6-pole, 1.5 m (501 16) for connecting the motion sensing element to the measurement transducer or to a CASSY-interface (e.g. 524 007) equipped with the BMW-box (524 032).

Appliance-plug connector on rear of device with integrated holder for primary and spare fuse. mains power lead and 2 fuses included in scope of supply (see section 4.2 for connection to mains voltage 115 V).

Two folding feet in bottom of housing for inclining the device.

Mains voltage: 230 V AC (set at factory), convertible to 115 V AC (see Section 4.2)  
 Power consumption: 6 VA  
 Dimensions: 20 cm x 15 cm x 15 cm  
 Weight: 0.93 kg

Scope of supply of motion transducer (337 63):

Motion sensing element (337 631) with measurement transducer  
 Multicore cable, 6-pole, 1.5 m (501 16)

### 3 Operation

#### 3.1 Setting up the motion sensing element

Use the stand material to set up the motion sensing element in the experiment arrangement so that the motion of the experiment body is transmitted with no slip via a thread to the spoked wheel (Fig. 2.1 - 2.3).

Lay the thread so that it runs in the plane of the thread groove.

#### 3.2 Measuring the voltage at the output of the measurement transducer (see Fig. 3.1/4)

In order to minimize measuring errors, use voltmeters with high accuracy ratings (class  $\leq 1.5$ ) and an internal resistance of  $\geq 20 \text{ k}\Omega$  at the s, v and a-outputs and  $\geq 200 \text{ k}\Omega$  at the t-output.

Recommended equipment:

AV-meter	531 94
or	
Demo-multimeter, zero point center	531 915
or	
Yt-recorder, 2-channel	575 712
or	
XY-Yt-recorder SR 720	575 663
or for computer-assisted processing of the output voltages of the measurement transducer	
CASSY-interface, e.g.	
CASSYpack-E	from 524 007
CASSY-Software	
- for DOS:	
Measuring and Evaluating	524 112
- for Windows:	
Universal Data Acquisition	525 032

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### 3.3 Synchronizing the start of the motion and the start of the time base at $t = 0$

Arrest the body whose motion is to be studied with a holding magnet (e.g. 336 21) which is connected to the interrupter sockets ⑬ via a voltage supply when switch ⑫ is in the reset position.

To start the motion and time base (time-proportional voltage at sockets ⑪) **quickly** set switch ⑫ to start; the connection across the interrupter socket pair ⑬ is interrupted when switch ⑫ is in the stop position, but the time base only starts at switch position "Start".

### 3.4 Operate delay; limits of measuring accuracy

Settling time of filters  $< 0.2$  s, i.e. for abrupt changes in the quantity to be measured the corresponding output voltage typically follows after  $0.1$  s.

Time base without delay

At velocities  $< 0.05$  m/s, the digital measuring principle becomes apparent through unevenness of the  $a(t)$ -signal, at  $v < 0.02$  m/s in the  $v(t)$ -signal as well. Due to the unbalance of the spoked wheel ( $< 0.02$  mm), the  $a(t)$ -signal additionally demonstrates a ripple, which is dependent on the velocity and can reach a maximum amplitude of  $0.02$  m/s<sup>2</sup> at  $v \approx 0.5$  m/s.

Solution:

- Connect an RC-element ( $10$  k $\Omega$ ,  $10$   $\mu$ F) between the respective output of the measuring electronics and the recorder input. This signal smoothing increases the settling time.
- Use a demonstration meter instead of the recorder (appropriate for quantities which do not change over time).


### 3.5 Processing the digital output signals of the motion sensing element in a computer (see Fig. 5)

Equipment required:

CASSY-interface, e.g. from	
CASSYpack-E	524 007
BMW-box	524 032
Multicore cable, 6-pole, 1.5 m	501 16
CASSY-software	
- for DOS:	
Motion (BMW)	524 702
- for Windows:	
Universal Data Acquisition	525 032

For detailed information with experiment examples, see the manual for the CASSY-software "Motion, BMW" (524 702) and the experiment examples contained on the program disks (524 702 and 525 032).

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Motion on the linear air track (337 50)

Display of the characteristic quantities path  $s$ , time  $t$ , velocity  $v$  and acceleration  $a$  of the motion of the trolley using a voltmeter connected to the outputs of the measurement transducer (from 337 63) (e.g.  $a = \text{const.}$  for uniformly accelerated motion)

Fig. 3.2

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Computer-assisted processing of signals from the motion sensing element (337 631) via a CASSY interface device (e.g. from 524 007) equipped with a BMW box (524 032): real-time registration of  $s(t)$  with e.g. display of  $s(t)$ ,  $v(t)$  and  $a(t)$  in a single diagram.


Fig. 4

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Recording the  $s(t)$ -,  $v(t)$  and  $a(t)$  diagrams of the oscillation of a helical spring with an XY-recorder using the time base of the measurement transducer (from 337 63)


Fig. 5

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Computer-assisted recording of harmonic and chaotic oscillations of the torsion balance (346 00)

Realizing the non-linear force field required for chaotic motion by asymmetrically attaching weights from 2 g to 20 g (e.g. from 590 27) to the pendulum body (e.g. from 530 27) to the pendulum body at the position A shown here using two-sided adhesive tape.

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
Fig. 6.1


Fig. 6.2

## 5 Replacing the fuse, converting the mains voltage

### 6 Replacing the primary fuse

Pry out insert (a) with holder for primary fuse (b) and spare fuse (c) (Fig. 6.1).

Replace the defective fuse (b) using a new fuse (c) which has been checked for the correct rating (Fig. 6.2).

Insert the new spare fuse (see at fuse label on rear of device for fuserating) (c) and put the insert (a) back in the device.

### 6.1 Converting the device for operation with 115 V mains voltage

*Attention: always unplug the device first!*

Loosen housing screw (e) on the bottom of the device using a Philips screwdriver (size 2) (Fig. 7.1).

Turn the device right side up and remove the top housing section (f) (Fig. 7.2).

Remove the PC board (h) in the 230 V position from the blue slot (g) at the transformer (Fig. 7.3).

Turn the board (h) around so that the stamped designation "110 V" (for 115 V mains voltage) is at the bottom left above the "1" on the contact side of the plug-in slot (g) (Fig. 7.4).

Insert board (h) and screw the housing back together.

Insert the correct primary fuse for the new mains voltage (see at fuse label on rear of device for fuserating).

See section 4.1 (Fig. 6.1 and 6.2) for instructions on how to change the fuse.


Fig. 7.1 - 7.4